

CWC15 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50720-4

Capture Antibody Information

Catalog Number:
60514-3-PBS
Host:
Mouse
Isotype:
IgG2b
Purification Method:
Protein A Magarose purification

Clone ID:
1A4E11
Reactivity:
human
GenBank:
BC040946
Immunogen Catalog Number:
Ag18323

Conjugate:
Unconjugated
Full name:
CWC15 spliceosome-associated protein homolog (S. cerevisiae)
Gene ID:
51503

Detection Antibody Information

Catalog Number:
60514-4-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
3C12G2
Reactivity:
human
GenBank:
BC040946
Immunogen Catalog Number:
Ag18323

Conjugate:
Unconjugated
Full name:
CWC15 spliceosome-associated protein homolog (S. cerevisiae)
Gene ID:
51503

Applications

Tested Applications:
Cytometric bead array

Range:
3.125-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.


Product Information

MP50720-4 targets CWC15 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: CWC15 Monoclonal antibody, PBS Only (Capture) 60514-3-PBS (1A4E11). 100 µg. Concentration 1 mg/ml.
Detection antibody: CWC15 Monoclonal antibody, PBS Only (Detector) 60514-4-PBS (3C12G2). 100 µg. Concentration 1 mg/ml.
Alternative CWC15 matched antibody pairs: MP50720-1, MP50720-2, MP50720-3
Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.
Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50720-4, CWC15 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60514-3-PBS. Detection antibody: 60514-4-PBS. Standard: Ag18323. Range: 3.125-100 ng/mL.