

# LACTB2 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50726-4

## Capture Antibody Information

<b>Catalog Number:</b> 67399-1-PBS	<b>Clone ID:</b> 1G6B1	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Mouse	<b>Reactivity:</b> human, mouse	<b>Full name:</b> lactamase, beta 2
<b>Isotype:</b> IgG1	<b>GenBank:</b> BC000878	<b>Gene ID:</b> 51110
<b>Purification Method:</b> Protein G purification	<b>Immunogen Catalog Number:</b> Ag10231	

## Detection Antibody Information

<b>Catalog Number:</b> 67399-4-PBS	<b>Clone ID:</b> 1C2G11	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Mouse	<b>Reactivity:</b> human	<b>Full name:</b> lactamase, beta 2
<b>Isotype:</b> IgG1	<b>GenBank:</b> BC000878	<b>Gene ID:</b> 51110
<b>Purification Method:</b> Protein G Magarose purification	<b>Immunogen Catalog Number:</b> Ag10231	

## Applications

<b>Tested Applications:</b> Cytometric bead array	<b>Range:</b> 3.125-100 ng/mL (Cytometric Bead Array)	<b>Recommended Dilutions:</b> It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

## Product Information

MP50726-4 targets LACTB2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: LACTB2 Monoclonal antibody, PBS Only (Capture) 67399-1-PBS (1G6B1). 100 µg. Concentration 1 mg/ml.

Detection antibody: LACTB2 Monoclonal antibody, PBS Only (Detector) 67399-4-PBS (1C2G11). 100 µg. Concentration 1 mg/ml.

Alternative LACTB2 matched antibody pairs: MP50726-1, MP50726-2, MP50726-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP50726-4, LACTB2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67399-1-PBS. Detection antibody: 67399-4-PBS. Standard: Ag10231. Range: 3.125-100 ng/mL.