

For Research Use Only

BCL7B Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50729-4

Capture Antibody Information

Catalog Number: 60517-5-PBS	Clone ID: 1D1H8	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: B-cell CLL/lymphoma 7B
Isotype: IgG1	GenBank: BC000956	Gene ID: 9275
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag22892	

Detection Antibody Information

Catalog Number: 60517-2-PBS	Clone ID: 3B11F11	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: B-cell CLL/lymphoma 7B
Isotype: IgG2b	GenBank: BC000956	Gene ID: 9275
Purification Method: Protein A Magarose purification	Immunogen Catalog Number: Ag22892	

Applications

Tested Applications: Cytometric bead array	Range: 0.195-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP50729-4 targets BCL7B in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: BCL7B Monoclonal antibody, PBS Only (Capture) 60517-5-PBS (1D1H8). 100 µg. Concentration 1 mg/ml.

Detection antibody: BCL7B Monoclonal antibody, PBS Only (Detector) 60517-2-PBS (3B11F11). 100 µg. Concentration 1 mg/ml.

Alternative BCL7B matched antibody pairs: MP50729-1, MP50729-2, MP50729-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50729-4, BCL7B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60517-5-PBS. Detection antibody: 60517-2-PBS. Standard: Ag22892. Range: 0.195-100 ng/mL.