

For Research Use Only

AEG-1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50732-4

Capture Antibody Information

Catalog Number:
68592-6-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1F5E2

Reactivity:
human

GenBank:
BC045642

Immunogen Catalog Number:
Ag5033

Conjugate:
Unconjugated

Full name:
metadherin

Gene ID:
92140

Detection Antibody Information

Catalog Number:
68592-4-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
2B11F8

Reactivity:
human

GenBank:
BC045642

Immunogen Catalog Number:
Ag5033

Conjugate:
Unconjugated

Full name:
metadherin

Gene ID:
92140

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50732-4 targets AEG-1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: AEG-1 Monoclonal antibody, PBS Only (Capture) 68592-6-PBS (1F5E2). 100 μ g. Concentration 1 mg/ml.

Detection antibody: AEG-1 Monoclonal antibody, PBS Only (Detector) 68592-4-PBS (2B11F8). 100 μ g. Concentration 1 mg/ml.

Alternative AEG-1 matched antibody pairs: MP50732-1, MP50732-2, MP50732-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50732-4, AEG-1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68592-6-PBS. Detection antibody: 68592-4-PBS. Standard: Ag5033. Range: 0.781-100 ng/mL.