

For Research Use Only

SGPL1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50739-1

Capture Antibody Information

Catalog Number:
60523-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
1D5B8
Reactivity:
human
GenBank:
NM_003901
Immunogen Catalog Number:
Ag26959

Conjugate:
Unconjugated
Full name:
sphingosine-1-phosphate lyase 1
Gene ID:
8879

Detection Antibody Information

Catalog Number:
60523-2-PBS
Host:
Mouse
Isotype:
IgG2a
Purification Method:
Protein A purification

Clone ID:
4F12C12
Reactivity:
human
GenBank:
NM_003901
Immunogen Catalog Number:
Ag26959

Conjugate:
Unconjugated
Full name:
sphingosine-1-phosphate lyase 1
Gene ID:
8879

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50739-1 targets SGPL1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: SGPL1 Monoclonal antibody, PBS Only (Capture) 60523-1-PBS (1D5B8). 100 µg. Concentration 1 mg/ml.
Detection antibody: SGPL1 Monoclonal antibody, PBS Only (Detector) 60523-2-PBS (4F12C12). 100 µg. Concentration 1 mg/ml.
Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.
Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50739-1, SGPL1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60523-1-PBS. Detection antibody: 60523-2-PBS. Standard: Ag26959. Range: 1.563-200 ng/mL.