

For Research Use Only

Calnexin Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50746-1

Capture Antibody Information

Catalog Number:
66903-2-PBS

Host:
Mouse

Isotype:
IgG2a

Purification Method:
Protein A Magarose purification

Clone ID:
1G11C1

Reactivity:
human

GenBank:
BC003552

Immunogen Catalog Number:
Ag25535

Conjugate:
Unconjugated

Full name:
calnexin

Gene ID:
821

Detection Antibody Information

Catalog Number:
66903-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1C6B6

Reactivity:
human

GenBank:
BC003552

Immunogen Catalog Number:
Ag25535

Conjugate:
Unconjugated

Full name:
calnexin

Gene ID:
821

Applications

Tested Applications:
Cytometric bead array

Range:
0.195-50 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50746-1 targets Calnexin in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Calnexin Monoclonal antibody, PBS Only (Capture) 66903-2-PBS (1G11C1). 100 µg. Concentration 1 mg/ml.

Detection antibody: Calnexin Monoclonal antibody, PBS Only (Detector) 66903-3-PBS (1C6B6). 100 µg. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50746-1, Calnexin Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66903-2-PBS. Detection antibody: 66903-3-PBS. Standard: Ag25535. Range: 0.195-50 ng/mL