

# COX6B2 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50751-3

## Capture Antibody Information

<b>Catalog Number:</b> 60531-1-PBS	<b>Clone ID:</b> 1A3B5	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Mouse	<b>Reactivity:</b> human	<b>Full name:</b> cytochrome c oxidase subunit VIb polypeptide 2 (testis)
<b>Isotype:</b> IgG1	<b>GenBank:</b> BC064548	<b>Gene ID:</b> 125965
<b>Purification Method:</b> Protein G Magarose purification	<b>Immunogen Catalog Number:</b> Ag33133	

## Detection Antibody Information

<b>Catalog Number:</b> 60531-4-PBS	<b>Clone ID:</b> 1E12C1	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Mouse	<b>Reactivity:</b> human	<b>Full name:</b> cytochrome c oxidase subunit VIb polypeptide 2 (testis)
<b>Isotype:</b> IgG1	<b>GenBank:</b> BC064548	<b>Gene ID:</b> 125965
<b>Purification Method:</b> Protein G Magarose purification	<b>Immunogen Catalog Number:</b> Ag33133	

## Applications

<b>Tested Applications:</b> Cytometric bead array	<b>Range:</b> 0.781-100 ng/mL (Cytometric Bead Array)	<b>Recommended Dilutions:</b> It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

## Product Information

MP50751-3 targets COX6B2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: COX6B2 Monoclonal antibody, PBS Only (Capture) 60531-1-PBS (1A3B5). 100 µg. Concentration 1 mg/ml.

Detection antibody: COX6B2 Monoclonal antibody, PBS Only (Detector) 60531-4-PBS (1E12C1). 100 µg. Concentration 1 mg/ml.

Alternative COX6B2 matched antibody pairs: MP00220-1, MP50751-1, MP50751-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP50751-3, COX6B2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60531-1-PBS. Detection antibody: 60531-4-PBS. Standard: Ag33133. Range: 0.781-100 ng/mL.