

# GREB1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50752-2

## Capture Antibody Information

**Catalog Number:**  
60532-3-PBS

**Host:**  
Mouse

**Isotype:**  
IgG2a

**Purification Method:**  
Protein A Magarose purification

**Clone ID:**  
1C1E4

**Reactivity:**  
human

**GenBank:**  
NM\_014668

**Immunogen Catalog Number:**  
Ag29941

**Conjugate:**  
Unconjugated

**Full name:**  
GREB1 protein

**Gene ID:**  
9687

## Detection Antibody Information

**Catalog Number:**  
60532-1-PBS

**Host:**  
Mouse

**Isotype:**  
IgG2b

**Purification Method:**  
Protein A purification

**Clone ID:**  
2F3C7

**Reactivity:**  
human

**GenBank:**  
NM\_014668

**Immunogen Catalog Number:**  
Ag29941

**Conjugate:**  
Unconjugated

**Full name:**  
GREB1 protein

**Gene ID:**  
9687

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.098-100 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50752-2 targets GREB1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: GREB1 Monoclonal antibody, PBS Only (Capture) 60532-3-PBS (1C1E4). 100 µg. Concentration 1 mg/mL.

Detection antibody: GREB1 Monoclonal antibody, PBS Only (Capture/Detector) 60532-1-PBS (2F3C7). 100 µg. Concentration 1 mg/mL.

Alternative GREB1 matched antibody pairs: MP50752-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP50752-2, GREB1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60532-3-PBS. Detection antibody: 60532-1-PBS. Standard: Ag29941. Range: 0.098-100 ng/mL.