

For Research Use Only

# DMWD Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50758-1

## Capture Antibody Information

**Catalog Number:**  
60536-1-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G Magarose purification

**Clone ID:**  
2G10B9  
**Reactivity:**  
human  
**GenBank:**  
BC019266  
**Immunogen Catalog Number:**  
Ag36845

**Conjugate:**  
Unconjugated  
**Full name:**  
dystrophia myotonica, WD repeat containing  
**Gene ID:**  
1762

## Detection Antibody Information

**Catalog Number:**  
60536-2-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G Magarose purification

**Clone ID:**  
1F2A2  
**Reactivity:**  
human  
**GenBank:**  
BC019266  
**Immunogen Catalog Number:**  
Ag36845

**Conjugate:**  
Unconjugated  
**Full name:**  
dystrophia myotonica, WD repeat containing  
**Gene ID:**  
1762

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.098-100 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50758-1 targets DMWD in immunoassays as a matched antibody pair. Validated in Cytometric bead array.  
Capture antibody: DMWD Monoclonal antibody, PBS Only (Capture) 60536-1-PBS (2G10B9). 100 µg. Concentration 1 mg/ml.  
Detection antibody: DMWD Monoclonal antibody, PBS Only (Capture/Detector) 60536-2-PBS (1F2A2). 100 µg. Concentration 1 mg/ml.  
Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.  
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.  
Antibody use should be optimized for each application and assay.


## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:  
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50758-1, DMWD Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60536-1-PBS. Detection antibody: 60536-2-PBS. Standard: Ag36845. Range: 0.098-100 ng/mL.