

For Research Use Only

NTF3 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50760-3

Capture Antibody Information

Catalog Number:
60520-4-PBS
Host:
Mouse
Isotype:
IgG2a
Purification Method:
Protein A Magarose purification

Clone ID:
2C1H3
Reactivity:
human
GenBank:
BC069773
Immunogen Catalog Number:
Ag35836

Conjugate:
Unconjugated
Full name:
neurotrophin 3
Gene ID:
4908

Detection Antibody Information

Catalog Number:
60520-1-PBS
Host:
Mouse
Isotype:
IgG3
Purification Method:
Protein A Magarose purification

Clone ID:
3G12F11
Reactivity:
human
GenBank:
BC069773
Immunogen Catalog Number:
Ag35836

Conjugate:
Unconjugated
Full name:
neurotrophin 3
Gene ID:
4908

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50760-3 targets NTF3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: NTF3 Monoclonal antibody, PBS Only (Capture) 60520-4-PBS (2C1H3). 100 µg. Concentration 1 mg/ml.
Detection antibody: NTF3 Monoclonal antibody, PBS Only (Detector) 60520-1-PBS (3G12F11). 100 µg. Concentration 1 mg/ml.
Alternative NTF3 matched antibody pairs: MP50760-1, MP50760-2
Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.
Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50760-3, NTF3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60520-4-PBS. Detection antibody: 60520-1-PBS. Standard: Ag12790. Range: 1.563-200 ng/mL.