

For Research Use Only

RTN4IP1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50762-2

Capture Antibody Information

Catalog Number:
60538-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
1C4F11

Reactivity:
human

GenBank:
BC006399

Immunogen Catalog Number:
Ag8161

Conjugate:
Unconjugated

Full name:
reticulon 4 interacting protein 1

Gene ID:
84816

Detection Antibody Information

Catalog Number:
60538-4-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
2B11A3

Reactivity:
human

GenBank:
BC006399

Immunogen Catalog Number:
Ag8161

Conjugate:
Unconjugated

Full name:
reticulon 4 interacting protein 1

Gene ID:
84816

Applications

Tested Applications:
Cytometric bead array

Range:
0.195-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50762-2 targets RTN4IP1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RTN4IP1 Monoclonal antibody, PBS Only (Capture) 60538-3-PBS (1C4F11). 100 μ g. Concentration 1 mg/mL.

Detection antibody: RTN4IP1 Monoclonal antibody, PBS Only (Detector) 60538-4-PBS (2B11A3). 100 μ g. Concentration 1 mg/mL.

Alternative RTN4IP1 matched antibody pairs: MP50762-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50762-2, RTN4IP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60538-3-PBS. Detection antibody: 60538-4-PBS. Standard: Ag8161. Range: 0.195-100 ng/mL.