

SHC1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50763-2

Capture Antibody Information

Catalog Number:
67194-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1G5A11

Reactivity:
human

GenBank:
NM_183001

Immunogen Catalog Number:
Ag25173

Conjugate:
Unconjugated

Full name:
SHC (Src homology 2 domain containing) transforming protein 1

Gene ID:
6464

Detection Antibody Information

Catalog Number:
67194-4-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1H4G2

Reactivity:
human

GenBank:
NM_183001

Immunogen Catalog Number:
Ag25173

Conjugate:
Unconjugated

Full name:
SHC (Src homology 2 domain containing) transforming protein 1

Gene ID:
6464

Applications

Tested Applications:
Cytometric bead array

Range:
0.391-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50763-2 targets SHC1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SHC1 Monoclonal antibody, PBS Only (Capture/Detector) 67194-3-PBS (1G5A11). 100 μ g. Concentration 1 mg/ml.

Detection antibody: SHC1 Monoclonal antibody, PBS Only (Detector) 67194-4-PBS (1H4G2). 100 μ g. Concentration 1 mg/ml.

Alternative SHC1 matched antibody pairs: MP50763-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50763-2, SHC1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67194-3-PBS. Detection antibody: 67194-4-PBS. Standard: Ag25173. Range: 0.391-100 ng/mL.