

# SHC1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50763-4

## Capture Antibody Information

**Catalog Number:**  
67194-5-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG2a  
**Purification Method:**  
Protein A purification

**Clone ID:**  
2C11A4  
**Reactivity:**  
human  
**GenBank:**  
NM\_183001  
**Immunogen Catalog Number:**  
Ag25173

**Conjugate:**  
Unconjugated  
**Full name:**  
SHC (Src homology 2 domain containing) transforming protein 1  
**Gene ID:**  
6464

## Detection Antibody Information

**Catalog Number:**  
67194-3-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G Magarose purification

**Clone ID:**  
1G5A11  
**Reactivity:**  
human  
**GenBank:**  
NM\_183001  
**Immunogen Catalog Number:**  
Ag25173

**Conjugate:**  
Unconjugated  
**Full name:**  
SHC (Src homology 2 domain containing) transforming protein 1  
**Gene ID:**  
6464

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.781-100 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50763-4 targets SHC1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SHC1 Monoclonal antibody, PBS Only (Capture) 67194-5-PBS (2C11A4). 100 µg. Concentration 1 mg/ml.

Detection antibody: SHC1 Monoclonal antibody, PBS Only (Capture/Detector) 67194-3-PBS (1G5A11). 100 µg. Concentration 1 mg/ml.

Alternative SHC1 matched antibody pairs: MP50763-1, MP50763-2, MP50763-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP50763-4, SHC1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67194-5-PBS. Detection antibody: 67194-3-PBS. Standard: Ag25173. Range: 0.781-100 ng/mL.