

For Research Use Only

ACSS1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50764-2

Capture Antibody Information

Catalog Number:
68608-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
3A9E10
Reactivity:
human
GenBank:
BC039261
Immunogen Catalog Number:
Ag10853

Conjugate:
Unconjugated
Full name:
acyl-CoA synthetase short-chain family member 1
Gene ID:
84532

Detection Antibody Information

Catalog Number:
68608-4-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
4A6C10
Reactivity:
human
GenBank:
BC039261
Immunogen Catalog Number:
Ag10853

Conjugate:
Unconjugated
Full name:
acyl-CoA synthetase short-chain family member 1
Gene ID:
84532

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50764-2 targets ACSS1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: ACSS1 Monoclonal antibody, PBS Only (Capture) 68608-2-PBS (3A9E10). 100 µg. Concentration 1 mg/mL.
Detection antibody: ACSS1 Monoclonal antibody, PBS Only (Detector) 68608-4-PBS (4A6C10). 100 µg. Concentration 1 mg/mL.
Alternative ACSS1 matched antibody pairs: MP50764-1
Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.
Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50764-2, ACS1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68608-2-PBS. Detection antibody: 68608-4-PBS. Standard: Ag10853. Range: 1.563-100 ng/mL.