

For Research Use Only

# STAR Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50772-2

## Capture Antibody Information

**Catalog Number:**  
67130-3-PBS

**Host:**  
Mouse

**Isotype:**  
IgG1

**Purification Method:**  
Protein G Magarose purification

**Clone ID:**  
1D5B1

**Reactivity:**  
human

**GenBank:**  
BC010550

**Immunogen Catalog Number:**  
Ag28698

**Conjugate:**  
Unconjugated

**Full name:**  
steroidogenic acute regulatory protein

**Gene ID:**  
6770

## Detection Antibody Information

**Catalog Number:**  
67130-4-PBS

**Host:**  
Mouse

**Isotype:**  
IgG1

**Purification Method:**  
Protein G Magarose purification

**Clone ID:**  
2E11A3

**Reactivity:**  
human

**GenBank:**  
BC010550

**Immunogen Catalog Number:**  
Ag28698

**Conjugate:**  
Unconjugated

**Full name:**  
steroidogenic acute regulatory protein

**Gene ID:**  
6770

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.195-200 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50772-2 targets STAR in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: STAR Monoclonal antibody, PBS Only (Capture/Detector) 67130-3-PBS (1D5B1). 100 µg. Concentration 1 mg/ml.

Detection antibody: STAR Monoclonal antibody, PBS Only (Detector) 67130-4-PBS (2E11A3). 100 µg. Concentration 1 mg/ml.

Alternative STAR matched antibody pairs: MP50772-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50772-2, STAR Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67130-3-PBS. Detection antibody: 67130-4-PBS. Standard: Ag28698. Range: 0.195-200 ng/mL.