

For Research Use Only

VAPA Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50773-2

Capture Antibody Information

Catalog Number:
67845-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1D2D8

Reactivity:
human

GenBank:
BC002992

Immunogen Catalog Number:
Ag7819

Conjugate:
Unconjugated

Full name:
VAMP (vesicle-associated membrane protein)-associated protein A, 33kDa

Gene ID:
9218

Detection Antibody Information

Catalog Number:
67845-4-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
1G12D7

Reactivity:
human

GenBank:
BC002992

Immunogen Catalog Number:
Ag7819

Conjugate:
Unconjugated

Full name:
VAMP (vesicle-associated membrane protein)-associated protein A, 33kDa

Gene ID:
9218

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50773-2 targets VAPA in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: VAPA Monoclonal antibody, PBS Only (Capture) 67845-2-PBS (1D2D8). 100 μ g. Concentration 1 mg/ml.

Detection antibody: VAPA Monoclonal antibody, PBS Only (Detector) 67845-4-PBS (1G12D7). 100 μ g. Concentration 1 mg/ml.

Alternative VAPA matched antibody pairs: MP50773-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50773-2, VAPA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67845-2-PBS. Detection antibody: 67845-4-PBS. Standard:Ag7819. Range: 1.563-200 ng/mL.