

# TET2 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50789-2

## Capture Antibody Information

<b>Catalog Number:</b> 60560-3-PBS	<b>Clone ID:</b> 1A1A12	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Mouse	<b>Reactivity:</b> human	<b>Full name:</b> tet oncogene family member 2
<b>Isotype:</b> IgG2b	<b>GenBank:</b> BC150180	<b>Gene ID:</b> 54790
<b>Purification Method:</b> Protein A Magarose purification	<b>Immunogen Catalog Number:</b> Ag15773	

## Detection Antibody Information

<b>Catalog Number:</b> 60560-4-PBS	<b>Clone ID:</b> 6B10C6	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Mouse	<b>Reactivity:</b> human	<b>Full name:</b> tet oncogene family member 2
<b>Isotype:</b> IgG1	<b>GenBank:</b> BC150180	<b>Gene ID:</b> 54790
<b>Purification Method:</b> Protein G purification	<b>Immunogen Catalog Number:</b> Ag15773	

## Applications

<b>Tested Applications:</b> Cytometric bead array	<b>Range:</b> 0.391-50 ng/mL (Cytometric Bead Array)	<b>Recommended Dilutions:</b> It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

## Product Information

MP50789-2 targets TET2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TET2 Monoclonal antibody, PBS Only (Capture) 60560-3-PBS (1A1A12). 100 µg. Concentration 1 mg/mL.

Detection antibody: TET2 Monoclonal antibody, PBS Only (Detector) 60560-4-PBS (6B10C6). 100 µg. Concentration 1 mg/mL.

Alternative TET2 matched antibody pairs: MP50789-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP50789-2, TET2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60560-3-PBS. Detection antibody: 60560-4-PBS. Standard: Ag15773. Range: 0.391-50 ng/mL