

SAR1B Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50803-3

Capture Antibody Information

Catalog Number: 67765-6-PBS	Clone ID: 1E3H5	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: SAR1 homolog B (<i>S. cerevisiae</i>)
Isotype: IgG2b	GenBank: BC002847	Gene ID: 51128
Purification Method: Protein A Magarose purification	Immunogen Catalog Number: Ag7718	

Detection Antibody Information

Catalog Number: 67765-5-PBS	Clone ID: 3B11C3	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: SAR1 homolog B (<i>S. cerevisiae</i>)
Isotype: IgG1	GenBank: BC002847	Gene ID: 51128
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag7718	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
------------------------------------------------------	----------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------

Product Information

MP50803-3 targets SAR1B in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SAR1B Monoclonal antibody, PBS Only (Capture) 67765-6-PBS (1E3H5). 100 μ g. Concentration 1 mg/mL.

Detection antibody: SAR1B Monoclonal antibody, PBS Only (Detector) 67765-5-PBS (3B11C3). 100 μ g. Concentration 1 mg/mL.

Alternative SAR1B matched antibody pairs: MP50803-1, MP50803-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50803-3, SAR1B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67765-6-PBS. Detection antibody: 67765-5-PBS. Standard: Ag7718. Range: 0.781-100 ng/mL.