

For Research Use Only

RENIN-RECEPTOR,ATP6AP2 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50807-1

Capture Antibody Information

Catalog Number:
60570-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
2B8G2
Reactivity:
human
GenBank:
BC010395
Immunogen Catalog Number:
Ag17725

Conjugate:
Unconjugated
Full name:
ATPase, H⁺ transporting, lysosomal accessory protein 2
Gene ID:
10159

Detection Antibody Information

Catalog Number:
60570-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
1F5E9
Reactivity:
human
GenBank:
BC010395
Immunogen Catalog Number:
Ag17725

Conjugate:
Unconjugated
Full name:
ATPase, H⁺ transporting, lysosomal accessory protein 2
Gene ID:
10159

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-25 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50807-1 targets RENIN-RECEPTOR,ATP6AP2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RENIN-RECEPTOR,ATP6AP2 Monoclonal antibody, PBS Only (Capture) 60570-1-PBS (2B8G2). 100 µg. Concentration 1 mg/ml.

Detection antibody: RENIN-RECEPTOR,ATP6AP2 Monoclonal antibody, PBS Only (Detector) 60570-2-PBS (1F5E9). 100 µg. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50807-1, RENIN-RECEPTOR,ATP6AP2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60570-1-PBS. Detection antibody: 60570-2-PBS. Standard:Ag17725. Range: 1.563-25 ng/mL