

For Research Use Only

PDIA6 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50818-4

Capture Antibody Information

Catalog Number:
66669-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1F12B3

Reactivity:
human

GenBank:
BC001312

Immunogen Catalog Number:
Ag27608

Conjugate:
Unconjugated

Full name:
protein disulfide isomerase family A,
member 6

Gene ID:
10130

Detection Antibody Information

Catalog Number:
66669-7-PBS

Host:
Mouse

Isotype:
IgG2a

Purification Method:
Protein A Magarose purification

Clone ID:
2A4E6

Reactivity:
human

GenBank:
BC001312

Immunogen Catalog Number:
Ag27608

Conjugate:
Unconjugated

Full name:
protein disulfide isomerase family A,
member 6

Gene ID:
10130

Applications

Tested Applications:
Cytometric bead array

Range:
3.125-100 ng/mL (Cytometric Bead
Array)

Recommended Dilutions:
It is recommended that this reagent
should be titrated in each testing
system to obtain optimal results.

Product Information

MP50818-4 targets PDIA6 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PDIA6 Monoclonal antibody, PBS Only (Capture) 66669-2-PBS (1F12B3). 100 μ g. Concentration 1 mg/mL.

Detection antibody: PDIA6 Monoclonal antibody, PBS Only (Detector) 66669-7-PBS (2A4E6). 100 μ g. Concentration 1 mg/mL.

Alternative PDIA6 matched antibody pairs: MP50818-1, MP50818-2, MP50818-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50818-4, PDIA6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66669-2-PBS. Detection antibody: 66669-7-PBS. Standard: Ag27608. Range: 3.125-100 ng/mL.