

For Research Use Only

IL15RA Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50819-1

Capture Antibody Information

Catalog Number:
60577-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
3C4G6
Reactivity:
human
GenBank:

Conjugate:
Unconjugated
Full name:
interleukin 15 receptor, alpha
Gene ID:
3601

Detection Antibody Information

Catalog Number:
60577-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
3A1C11
Reactivity:
human
GenBank:

Conjugate:
Unconjugated
Full name:
interleukin 15 receptor, alpha
Gene ID:
3601

Applications

Tested Applications:
Cytometric bead array

Range:
0.391-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50819-1 targets IL15RA in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: IL15RA Monoclonal antibody, PBS Only (Capture) 60577-1-PBS (3C4G6). 100 μ g. Concentration 1 mg/mL.
Detection antibody: IL15RA Monoclonal antibody, PBS Only (Detector) 60577-2-PBS (3A1C11). 100 μ g. Concentration 1 mg/mL.
Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.
Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50819-1, IL15RA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60577-1-PBS. Detection antibody: 60577-2-PBS. Standard:HZ-1323RA. Range: 0.391-100 ng/mL