

ADCY3 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: MP50823-2

Capture Antibody Information

Catalog Number: 60581-3-PBS	Clone ID: 1D11H3	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: adenylate cyclase 3
Isotype: IgG1	GenBank: BC126235	Gene ID: 109
Purification Method: Protein G purification	Immunogen Catalog Number: Ag25859	

Detection Antibody Information

Catalog Number: 60581-2-PBS	Clone ID: 2G1C1	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: adenylate cyclase 3
Isotype: IgG1	GenBank: BC126235	Gene ID: 109
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag25859	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP50823-2 targets ADCY3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ADCY3 Monoclonal antibody, PBS Only (Capture) 60581-3-PBS (1D11H3). 100 µg. Concentration 1 mg/mL.

Detection antibody: ADCY3 Monoclonal antibody, PBS Only (Detector) 60581-2-PBS (2G1C1). 100 µg. Concentration 1 mg/mL.

Alternative ADCY3 matched antibody pairs: MP50823-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50823-2, ADCY3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60581-3-PBS. Detection antibody: 60581-2-PBS. Standard: Ag25859. Range: 0.781-100 ng/mL.