

CREB3 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50828-1

Capture Antibody Information

Catalog Number:
68379-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
2F4H8
Reactivity:
human
GenBank:
BC009402
Immunogen Catalog Number:
Ag1823

Conjugate:
Unconjugated
Full name:
cAMP responsive element binding protein 3
Gene ID:
10488

Detection Antibody Information

Catalog Number:
68379-3-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
2C5E10
Reactivity:
human
GenBank:
BC009402
Immunogen Catalog Number:
Ag1823

Conjugate:
Unconjugated
Full name:
cAMP responsive element binding protein 3
Gene ID:
10488

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-25 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50828-1 targets CREB3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CREB3 Monoclonal antibody, PBS Only (Capture) 68379-2-PBS (2F4H8). 100 µg. Concentration 1 mg/ml.

Detection antibody: CREB3 Monoclonal antibody, PBS Only (Detector) 68379-3-PBS (2C5E10). 100 µg. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50828-1, CREB3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68379-2-PBS. Detection antibody: 68379-3-PBS. Standard: Ag1823. Range: 0.098-25 ng/mL