

For Research Use Only

# NFS1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50834-1

## Capture Antibody Information

**Catalog Number:**  
67021-2-PBS

**Host:**  
Mouse

**Isotype:**  
IgG1

**Purification Method:**  
Protein G Magarose purification

**Clone ID:**  
1F9A3

**Reactivity:**  
human

**GenBank:**  
BC018471

**Immunogen Catalog Number:**  
Ag8017

**Conjugate:**  
Unconjugated

**Full name:**  
NFS1 nitrogen fixation 1 homolog (S. cerevisiae)

**Gene ID:**  
9054

## Detection Antibody Information

**Catalog Number:**  
67021-3-PBS

**Host:**  
Mouse

**Isotype:**  
IgG1

**Purification Method:**  
Protein G Magarose purification

**Clone ID:**  
3A5E11

**Reactivity:**  
human

**GenBank:**  
BC018471

**Immunogen Catalog Number:**  
Ag8017

**Conjugate:**  
Unconjugated

**Full name:**  
NFS1 nitrogen fixation 1 homolog (S. cerevisiae)

**Gene ID:**  
9054

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.781-200 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50834-1 targets NFS1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NFS1 Monoclonal antibody, PBS Only (Capture) 67021-2-PBS (1F9A3). 100  $\mu$ g. Concentration 1 mg/ml.

Detection antibody: NFS1 Monoclonal antibody, PBS Only (Detector) 67021-3-PBS (3A5E11). 100  $\mu$ g. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage


**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:  
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50834-1, NFS1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67021-2-PBS. Detection antibody: 67021-3-PBS. Standard: Ag8017. Range: 0.781-200 ng/mL.