

PTGIS Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50835-1

Capture Antibody Information

Catalog Number:
60586-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
1F10A9
Reactivity:
human
GenBank:
BC101809
Immunogen Catalog Number:
Ag25824

Conjugate:
Unconjugated
Full name:
prostaglandin I2 (prostacyclin) synthase
Gene ID:
5740

Detection Antibody Information

Catalog Number:
60586-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
1C7F2
Reactivity:
human
GenBank:
BC101809
Immunogen Catalog Number:
Ag25824

Conjugate:
Unconjugated
Full name:
prostaglandin I2 (prostacyclin) synthase
Gene ID:
5740

Applications

Tested Applications:
Cytometric bead array

Range:
0.391-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50835-1 targets PTGIS in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PTGIS Monoclonal antibody, PBS Only (Capture) 60586-1-PBS (1F10A9). 100 µg. Concentration 1 mg/ml.

Detection antibody: PTGIS Monoclonal antibody, PBS Only (Detector) 60586-2-PBS (1C7F2). 100 µg. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50835-1, PTGIS Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60586-1-PBS. Detection antibody: 60586-2-PBS. Standard: Ag25824. Range: 0.391-100 ng/mL.