

For Research Use Only

RPS16 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50839-3

Capture Antibody Information

Catalog Number: 60590-1-PBS	Clone ID: 1B11H3	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: ribosomal protein S16
Isotype: IgG1	GenBank: BC004324	Gene ID: 6217
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag8385	

Detection Antibody Information

Catalog Number: 60590-4-PBS	Clone ID: 2D10F5	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: ribosomal protein S16
Isotype: IgG1	GenBank: BC004324	Gene ID: 6217
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag8385	

Applications

Tested Applications: Cytometric bead array	Range: 0.391-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP50839-3 targets RPS16 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RPS16 Monoclonal antibody, PBS Only (Capture) 60590-1-PBS (1B11H3). 100 µg. Concentration 1 mg/ml.

Detection antibody: RPS16 Monoclonal antibody, PBS Only (Detector) 60590-4-PBS (2D10F5). 100 µg. Concentration 1 mg/ml.

Alternative RPS16 matched antibody pairs: MP50839-1, MP50839-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50839-3, RPS16 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60590-1-PBS. Detection antibody: 60590-4-PBS. Standard: Ag8385. Range: 0.391-100 ng/mL.