

For Research Use Only

SAMD8 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50842-2

Capture Antibody Information

Catalog Number:
60592-1-PBS
Host:
Mouse
Isotype:
IgG2b
Purification Method:
Protein A purification

Clone ID:
1C7C7
Reactivity:
human
GenBank:
BC146556
Immunogen Catalog Number:
Ag33150

Conjugate:
Unconjugated
Full name:
sterile alpha motif domain containing 8
Gene ID:
142891

Detection Antibody Information

Catalog Number:
60592-3-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
2F1E11
Reactivity:
human
GenBank:
BC146556
Immunogen Catalog Number:
Ag33150

Conjugate:
Unconjugated
Full name:
sterile alpha motif domain containing 8
Gene ID:
142891

Applications

Tested Applications:
Cytometric bead array

Range:
6.25-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50842-2 targets SAMD8 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: SAMD8 Monoclonal antibody, PBS Only (Capture) 60592-1-PBS (1C7C7). 100 µg. Concentration 1 mg/mL.
Detection antibody: SAMD8 Monoclonal antibody, PBS Only (Detector) 60592-3-PBS (2F1E11). 100 µg. Concentration 1 mg/mL.
Alternative SAMD8 matched antibody pairs: MP50842-1
Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.
Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50842-2, SAMD8 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60592-1-PBS. Detection antibody: 60592-3-PBS. Standard:null. Range: null.