

For Research Use Only

# CDO1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50867-3

## Capture Antibody Information

<b>Catalog Number:</b> 60610-1-PBS	<b>Clone ID:</b> 1D2E3	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Mouse	<b>Reactivity:</b> human	<b>Full name:</b> cysteine dioxygenase, type I
<b>Isotype:</b> IgG1	<b>GenBank:</b> BC024241	<b>Gene ID:</b> 1036
<b>Purification Method:</b> Protein G Magarose purification	<b>Immunogen Catalog Number:</b> Ag3294	

## Detection Antibody Information

<b>Catalog Number:</b> 60610-4-PBS	<b>Clone ID:</b> 2B9A6	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Mouse	<b>Reactivity:</b> human	<b>Full name:</b> cysteine dioxygenase, type I
<b>Isotype:</b> IgG1	<b>GenBank:</b> BC024241	<b>Gene ID:</b> 1036
<b>Purification Method:</b> Protein G Magarose purification	<b>Immunogen Catalog Number:</b> Ag3294	

## Applications

<b>Tested Applications:</b> Cytometric bead array	<b>Range:</b> 0.098-100 ng/mL (Cytometric Bead Array)	<b>Recommended Dilutions:</b> It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

## Product Information

MP50867-3 targets CDO1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CDO1 Monoclonal antibody, PBS Only (Capture/Detector) 60610-1-PBS (1D2E3). 100 µg. Concentration 1 mg/ml.

Detection antibody: CDO1 Monoclonal antibody, PBS Only (Detector) 60610-4-PBS (2B9A6). 100 µg. Concentration 1 mg/ml.

Alternative CDO1 matched antibody pairs: MP50867-1, MP50867-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage


**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:  
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50867-3, CDO1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60610-1-PBS. Detection antibody: 60610-4-PBS. Standard:Ag3294. Range: 0.098-100 ng/mL.