

For Research Use Only

SILV Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50869-2

Capture Antibody Information

Catalog Number: 60612-3-PBS	Clone ID: 1A6B2	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: silver homolog (mouse)
Isotype: IgG1	GenBank: BC001414	Gene ID: 6490
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag6726	

Detection Antibody Information

Catalog Number: 60612-2-PBS	Clone ID: 1C8D6	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: silver homolog (mouse)
Isotype: IgG3	GenBank: BC001414	Gene ID: 6490
Purification Method: Protein A Magarose purification	Immunogen Catalog Number: Ag6726	

Applications

Tested Applications: Cytometric bead array	Range: 0.098-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP50869-2 targets SILV in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SILV Monoclonal antibody, PBS Only (Capture) 60612-3-PBS (1A6B2). 100 µg. Concentration 1 mg/mL.

Detection antibody: SILV Monoclonal antibody, PBS Only (Detector) 60612-2-PBS (1C8D6). 100 µg. Concentration 1 mg/mL.

Alternative SILV matched antibody pairs: MP50869-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50869-2, SILV Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60612-3-PBS. Detection antibody: 60612-2-PBS. Standard: Ag6726. Range: 0.098-100 ng/mL.