

# SILV Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50869-3

## Capture Antibody Information

<b>Catalog Number:</b> 60612-4-PBS	<b>Clone ID:</b> 2A6H6	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Mouse	<b>Reactivity:</b> human	<b>Full name:</b> silver homolog (mouse)
<b>Isotype:</b> IgG1	<b>GenBank:</b> BC001414	<b>Gene ID:</b> 6490
<b>Purification Method:</b> Protein G Magarose purification	<b>Immunogen Catalog Number:</b> Ag6726	

## Detection Antibody Information

<b>Catalog Number:</b> 60612-5-PBS	<b>Clone ID:</b> 1C10F8	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Mouse	<b>Reactivity:</b> human	<b>Full name:</b> silver homolog (mouse)
<b>Isotype:</b> IgG1	<b>GenBank:</b> BC001414	<b>Gene ID:</b> 6490
<b>Purification Method:</b> Protein G Magarose purification	<b>Immunogen Catalog Number:</b> Ag6726	

## Applications

<b>Tested Applications:</b> Cytometric bead array	<b>Range:</b> 0.098-6.25 ng/mL (Cytometric Bead Array)	<b>Recommended Dilutions:</b> It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

## Product Information

MP50869-3 targets SILV in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SILV Monoclonal antibody, PBS Only (Capture) 60612-4-PBS (2A6H6). 100 µg. Concentration 1 mg/ml.

Detection antibody: SILV Monoclonal antibody, PBS Only (Detector) 60612-5-PBS (1C10F8). 100 µg. Concentration 1 mg/ml.

Alternative SILV matched antibody pairs: MP50869-1, MP50869-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP50869-3, SILV Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60612-4-PBS. Detection antibody: 60612-5-PBS. Standard:Ag6726. Range: 0.098-6.25 ng/mL