

TEX264 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50870-2

Capture Antibody Information

Catalog Number: 68507-2-PBS	Clone ID: 1B5B10	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: testis expressed 264
Isotype: IgG1	GenBank: BC008742	Gene ID: 51368
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag33710	

Detection Antibody Information

Catalog Number: 68507-3-PBS	Clone ID: 1D8H8	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: testis expressed 264
Isotype: IgG1	GenBank: BC008742	Gene ID: 51368
Purification Method: Protein G purification	Immunogen Catalog Number: Ag33710	

Applications

Tested Applications: Cytometric bead array	Range: 0.098-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP50870-2 targets TEX264 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TEX264 Monoclonal antibody, PBS Only (Capture/Detector) 68507-2-PBS (1B5B10). 100 µg. Concentration 1 mg/ml.

Detection antibody: TEX264 Monoclonal antibody, PBS Only (Detector) 68507-3-PBS (1D8H8). 100 µg. Concentration 1 mg/ml.

Alternative TEX264 matched antibody pairs: MP50870-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50870-2, TEX264 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68507-2-PBS. Detection antibody: 68507-3-PBS. Standard: Ag33710. Range: 0.098-100 ng/mL.