

For Research Use Only

# KIAA0776 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50874-2

## Capture Antibody Information

Catalog Number:  
60614-3-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G purification

Clone ID:  
1D3A6

Reactivity:  
human

GenBank:  
BC036379

Immunogen Catalog Number:  
Ag23404

Conjugate:  
Unconjugated

Full name:  
KIAA0776

Gene ID:  
23376

## Detection Antibody Information

Catalog Number:  
60614-4-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G Magarose purification

Clone ID:  
2C10F8

Reactivity:  
human

GenBank:  
BC036379

Immunogen Catalog Number:  
Ag23404

Conjugate:  
Unconjugated

Full name:  
KIAA0776

Gene ID:  
23376

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.391-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50874-2 targets KIAA0776 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: KIAA0776 Monoclonal antibody, PBS Only (Capture) 60614-3-PBS (1D3A6). 100  $\mu$ g. Concentration 1 mg/ml.

Detection antibody: KIAA0776 Monoclonal antibody, PBS Only (Detector) 60614-4-PBS (2C10F8). 100  $\mu$ g. Concentration 1 mg/ml.

Alternative KIAA0776 matched antibody pairs: MP00135-1, MP00135-2, MP00135-3, MP50874-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50874-2, KIAA0776 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60614-3-PBS. Detection antibody: 60614-4-PBS. Standard: Ag23404. Range: 0.391-200 ng/mL