

GRAMD1B Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50885-3

Capture Antibody Information

| | | |
|--|---|--|
| Catalog Number: 67551-5-PBS | Clone ID: 1B12C5 | Conjugate: Unconjugated |
| Host: Mouse | Reactivity: human | Full name: GRAM domain containing 1B |
| Isotype: IgG2a | GenBank: BC107480 | Gene ID: 57476 |
| Purification Method: Protein A Magarose purification | Immunogen Catalog Number: Ag20670 | |

Detection Antibody Information

| | | |
|---|---|--|
| Catalog Number: 67551-6-PBS | Clone ID: 1B2H1 | Conjugate: Unconjugated |
| Host: Mouse | Reactivity: human | Full name: GRAM domain containing 1B |
| Isotype: IgG1 | GenBank: BC107480 | Gene ID: 57476 |
| Purification Method: Protein G purification | Immunogen Catalog Number: Ag20670 | |

Applications

| | | |
|--|---|---|
| Tested Applications: Cytometric bead array | Range: 0.098-25 ng/mL (Cytometric Bead Array) | Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results. |
|--|---|---|

Product Information

MP50885-3 targets GRAMD1B in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: GRAMD1B Monoclonal antibody, PBS Only (Capture) 67551-5-PBS (1B12C5). 100 µg. Concentration 1 mg/ml.

Detection antibody: GRAMD1B Monoclonal antibody, PBS Only (Detector) 67551-6-PBS (1B2H1). 100 µg. Concentration 1 mg/ml.

Alternative GRAMD1B matched antibody pairs: MP50885-1, MP50885-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50885-3, GRAMD1B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67551-5-PBS. Detection antibody: 67551-6-PBS. Standard: Ag20670. Range: 0.098-25 ng/mL