

GRAMD1B Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50885-4

Capture Antibody Information

Catalog Number: 67551-7-PBS	Clone ID: 1H9H2	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: GRAM domain containing 1B
Isotype: IgG2a	GenBank: BC107480	Gene ID: 57476
Purification Method: Protein A Magarose purification	Immunogen Catalog Number: Ag20670	

Detection Antibody Information

Catalog Number: 67551-3-PBS	Clone ID: 1C2F2	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: GRAM domain containing 1B
Isotype: IgG1	GenBank: BC107480	Gene ID: 57476
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag20670	

Applications

Tested Applications: Cytometric bead array	Range: 1.563-50 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP50885-4 targets GRAMD1B in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: GRAMD1B Monoclonal antibody, PBS Only (Capture) 67551-7-PBS (1H9H2). 100 µg. Concentration 1 mg/ml.

Detection antibody: GRAMD1B Monoclonal antibody, PBS Only (Detector) 67551-3-PBS (1C2F2). 100 µg. Concentration 1 mg/ml.

Alternative GRAMD1B matched antibody pairs: MP50885-1, MP50885-2, MP50885-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50885-4, GRAMD1B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67551-7-PBS. Detection antibody: 67551-3-PBS. Standard: Ag20670. Range: 1.563-50 ng/mL