

For Research Use Only

RABL3 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50900-2

Capture Antibody Information

Catalog Number:
68342-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
2G2G2
Reactivity:
human
GenBank:
BC020832
Immunogen Catalog Number:
Ag33137

Conjugate:
Unconjugated
Full name:
RAB, member of RAS oncogene family-like 3
Gene ID:
285282

Detection Antibody Information

Catalog Number:
68342-4-PBS
Host:
Mouse
Isotype:
IgG2b
Purification Method:
Protein A Magarose purification

Clone ID:
2F4E10
Reactivity:
human
GenBank:
BC020832
Immunogen Catalog Number:
Ag33137

Conjugate:
Unconjugated
Full name:
RAB, member of RAS oncogene family-like 3
Gene ID:
285282

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50900-2 targets RABL3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: RABL3 Monoclonal antibody, PBS Only (Capture) 68342-2-PBS (2G2G2). 100 μ g. Concentration 1 mg/ml.
Detection antibody: RABL3 Monoclonal antibody, PBS Only (Detector) 68342-4-PBS (2F4E10). 100 μ g. Concentration 1 mg/ml.
Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.
Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50900-2, RABL3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68342-2-PBS. Detection antibody: 68342-4-PBS. Standard: Ag33137. Range: 0.098-100 ng/mL.