

For Research Use Only

FMO4 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50902-2

Capture Antibody Information

Catalog Number:
67881-4-PBS

Host:
Mouse

Isotype:
IgG3

Purification Method:
Protein A Magarose purification

Clone ID:
1B12B4

Reactivity:
human

GenBank:
BC002780

Immunogen Catalog Number:
Ag7935

Conjugate:
Unconjugated

Full name:
flavin containing monooxygenase 4

Gene ID:
2329

Detection Antibody Information

Catalog Number:
67881-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1B5H4

Reactivity:
human

GenBank:
BC002780

Immunogen Catalog Number:
Ag7935

Conjugate:
Unconjugated

Full name:
flavin containing monooxygenase 4

Gene ID:
2329

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50902-2 targets FMO4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FMO4 Monoclonal antibody, PBS Only (Capture) 67881-4-PBS (1B12B4). 100 μ g. Concentration 1 mg/ml.

Detection antibody: FMO4 Monoclonal antibody, PBS Only (Capture/Detector) 67881-2-PBS (1B5H4). 100 μ g. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50902-2, FMO4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67881-4-PBS. Detection antibody: 67881-2-PBS. Standard: Ag7935. Range: 0.098-100 ng/mL.