

For Research Use Only

MBTPS1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50910-2

Capture Antibody Information

Catalog Number:
60636-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
2G10A2

Reactivity:
human

GenBank:
BC026330

Immunogen Catalog Number:
Ag26396

Conjugate:
Unconjugated

Full name:
membrane-bound transcription factor
peptidase, site 1

Gene ID:
8720

Detection Antibody Information

Catalog Number:
60636-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
2A1B11

Reactivity:
human

GenBank:
BC026330

Immunogen Catalog Number:
Ag26396

Conjugate:
Unconjugated

Full name:
membrane-bound transcription factor
peptidase, site 1

Gene ID:
8720

Applications

Tested Applications:
Cytometric bead array

Range:
3.125-200 ng/mL (Cytometric Bead
Array)

Recommended Dilutions:
It is recommended that this reagent
should be titrated in each testing
system to obtain optimal results.

Product Information

MP50910-2 targets MBTPS1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MBTPS1 Monoclonal antibody, PBS Only (Capture) 60636-3-PBS (2G10A2). 100 μ g. Concentration 1 mg/ml.

Detection antibody: MBTPS1 Monoclonal antibody, PBS Only (Detector) 60636-2-PBS (2A1B11). 100 μ g. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50910-2, MBTPS1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60636-3-PBS. Detection antibody: 60636-2-PBS. Standard:Ag26396. Range: 3.125-200 ng/mL.