

For Research Use Only

GPR177 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50915-2

Capture Antibody Information

Catalog Number:
67435-4-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
1B6A6

Reactivity:
human

GenBank:
BC110826

Immunogen Catalog Number:
Ag12446

Conjugate:
Unconjugated

Full name:
G protein-coupled receptor 177

Gene ID:
79971

Detection Antibody Information

Catalog Number:
67435-5-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1D1H8

Reactivity:
human

GenBank:
BC110826

Immunogen Catalog Number:
Ag12446

Conjugate:
Unconjugated

Full name:
G protein-coupled receptor 177

Gene ID:
79971

Applications

Tested Applications:
Cytometric bead array

Range:
6.25-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50915-2 targets GPR177 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: GPR177 Monoclonal antibody, PBS Only (Capture) 67435-4-PBS (1B6A6). 100 μ g. Concentration 1 mg/ml.

Detection antibody: GPR177 Monoclonal antibody, PBS Only (Detector) 67435-5-PBS (1D1H8). 100 μ g. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50915-2, GPR177 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67435-4-PBS. Detection antibody: 67435-5-PBS. Standard: Ag12446. Range: 6.25-200 ng/mL