

For Research Use Only

EDN3 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50920-2

Capture Antibody Information

Catalog Number:
60643-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
1G11A5
Reactivity:
human
GenBank:
BC008876
Immunogen Catalog Number:
Ag1071

Conjugate:
Unconjugated
Full name:
endothelin 3
Gene ID:
1908

Detection Antibody Information

Catalog Number:
60642-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
1G12A11
Reactivity:
human
GenBank:
BC008876
Immunogen Catalog Number:
Ag26261

Conjugate:
Unconjugated
Full name:
endothelin 3
Gene ID:
1908

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50920-2 targets EDN3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: EDN3 Monoclonal antibody, PBS Only (Capture) 60643-2-PBS (1G11A5). 100 µg. Concentration 1 mg/ml.

Detection antibody: EDN3 Monoclonal antibody, PBS Only (Detector) 60642-2-PBS (1G12A11). 100 µg. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50920-2, EDN3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60643-2-PBS. Detection antibody: 60642-2-PBS. Standard: Ag1071. Range: 0.098-100 ng/mL.