

# KIAA0319L Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50926-1

## Capture Antibody Information

**Catalog Number:**  
68717-2-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G Magarose purification

**Clone ID:**  
2D10C6  
**Reactivity:**  
human  
**GenBank:**  
BC031672  
**Immunogen Catalog Number:**  
Ag15263

**Conjugate:**  
Unconjugated  
**Full name:**  
KIAA0319-like  
**Gene ID:**  
79932

## Detection Antibody Information

**Catalog Number:**  
68717-3-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG2b  
**Purification Method:**  
Protein A Magarose purification

**Clone ID:**  
2B10G12  
**Reactivity:**  
human  
**GenBank:**  
BC031672  
**Immunogen Catalog Number:**  
Ag15263

**Conjugate:**  
Unconjugated  
**Full name:**  
KIAA0319-like  
**Gene ID:**  
79932

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.098-100 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50926-1 targets KIAA0319L in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: KIAA0319L Monoclonal antibody, PBS Only (Capture) 68717-2-PBS (2D10C6). 100 µg. Concentration 1 mg/ml.

Detection antibody: KIAA0319L Monoclonal antibody, PBS Only (Detector) 68717-3-PBS (2B10G12). 100 µg. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


## Storage

**Storage:**  
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP50926-1, KIAA0319L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68717-2-PBS. Detection antibody: 68717-3-PBS. Standard: Ag15263. Range: 0.098-100 ng/mL