

For Research Use Only

# AIFM3 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50942-1

## Capture Antibody Information

**Catalog Number:**  
60658-1-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G purification

**Clone ID:**  
3C5G3  
**Reactivity:**  
human  
**GenBank:**  
BC032485  
**Immunogen Catalog Number:**  
Ag6278

**Conjugate:**  
Unconjugated  
**Full name:**  
apoptosis-inducing factor, mitochondrion-associated, 3  
**Gene ID:**  
150209

## Detection Antibody Information

**Catalog Number:**  
60658-2-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G purification

**Clone ID:**  
2B5D3  
**Reactivity:**  
human  
**GenBank:**  
BC032485  
**Immunogen Catalog Number:**  
Ag6278

**Conjugate:**  
Unconjugated  
**Full name:**  
apoptosis-inducing factor, mitochondrion-associated, 3  
**Gene ID:**  
150209

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.781-200 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50942-1 targets AIFM3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: AIFM3 Monoclonal antibody, PBS Only (Capture) 60658-1-PBS (3C5G3). 100 µg. Concentration 1 mg/ml.

Detection antibody: AIFM3 Monoclonal antibody, PBS Only (Detector) 60658-2-PBS (2B5D3). 100 µg. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50942-1, AIFM3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60658-1-PBS. Detection antibody: 60658-2-PBS. Standard:Ag6278. Range: 0.781-200 ng/mL.