

For Research Use Only

C3orf10 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50946-4

Capture Antibody Information

Catalog Number:
60663-4-PBS

Host:
Mouse

Isotype:
IgG2a

Purification Method:
Protein A Magarose purification

Clone ID:
2E9D7

Reactivity:
human

GenBank:
BC019303

Immunogen Catalog Number:
Ag33165

Conjugate:
Unconjugated

Full name:
chromosome 3 open reading frame 10

Gene ID:
55845

Detection Antibody Information

Catalog Number:
60663-7-PBS

Host:
Mouse

Isotype:
IgG2a

Purification Method:
Protein A Magarose purification

Clone ID:
1G3B10

Reactivity:
human

GenBank:
BC019303

Immunogen Catalog Number:
Ag33165

Conjugate:
Unconjugated

Full name:
chromosome 3 open reading frame 10

Gene ID:
55845

Applications

Tested Applications:
Cytometric bead array

Range:
0.195-50 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50946-4 targets C3orf10 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: C3orf10 Monoclonal antibody, PBS Only (Capture/Detector) 60663-4-PBS (2E9D7). 100 μ g. Concentration 1 mg/ml.

Detection antibody: C3orf10 Monoclonal antibody, PBS Only (Detector) 60663-7-PBS (1G3B10). 100 μ g. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50946-4, C3orf10 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60663-4-PBS. Detection antibody: 60663-7-PBS. Standard:Ag33165. Range: 0.195-50 ng/mL