

For Research Use Only

CSNK2A1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50947-2

Capture Antibody Information

Catalog Number:
68200-4-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
2B2D10

Reactivity:
human

GenBank:
BC011668

Immunogen Catalog Number:
Ag29062

Conjugate:
Unconjugated

Full name:
casein kinase 2, alpha 1 polypeptide

Gene ID:
1457

Detection Antibody Information

Catalog Number:
68200-5-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
3B2G7

Reactivity:
human

GenBank:
BC011668

Immunogen Catalog Number:
Ag29062

Conjugate:
Unconjugated

Full name:
casein kinase 2, alpha 1 polypeptide

Gene ID:
1457

Applications

Tested Applications:
Cytometric bead array

Range:
6.25-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50947-2 targets CSNK2A1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CSNK2A1 Monoclonal antibody, PBS Only (Capture) 68200-4-PBS (2B2D10). 100 μ g. Concentration 1 mg/ml.

Detection antibody: CSNK2A1 Monoclonal antibody, PBS Only (Detector) 68200-5-PBS (3B2G7). 100 μ g. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50947-2, CSNK2A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68200-4-PBS. Detection antibody: 68200-5-PBS. Standard: Ag29062. Range: 6.25-100 ng/mL.