

For Research Use Only

DNM2 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50954-1

Capture Antibody Information

Catalog Number:
68209-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein A purification

Clone ID:
1B8H3
Reactivity:
human, mouse, rat, pig, rabbit
GenBank:
BC054501
Immunogen Catalog Number:
Ag6274

Conjugate:
Unconjugated
Full name:
dynamin 2
Gene ID:
1785

Detection Antibody Information

Catalog Number:
68209-2-PBS
Host:
Mouse
Isotype:
IgG3
Purification Method:
Protein A Magarose purification

Clone ID:
1A5F9
Reactivity:
human
GenBank:
BC054501
Immunogen Catalog Number:
Ag6274

Conjugate:
Unconjugated
Full name:
dynamin 2
Gene ID:
1785

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50954-1 targets DNM2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: DNM2 Monoclonal antibody, PBS Only (Capture) 68209-1-PBS (1B8H3). 100 µg. Concentration 1 mg/ml.

Detection antibody: DNM2 Monoclonal antibody, PBS Only (Detector) 68209-2-PBS (1A5F9). 100 µg. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50954-1, DNM2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68209-1-PBS. Detection antibody: 68209-2-PBS. Standard:Ag6274. Range: 0.781-100 ng/mL.