

For Research Use Only

# OTUD6B Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50962-2

## Capture Antibody Information

Catalog Number:  
68450-1-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein A purification

Clone ID:  
1D1G1

Reactivity:  
human, mouse, rat

GenBank:  
BC029760

Immunogen Catalog Number:  
Ag33607

Conjugate:

Unconjugated

Full name:

OTU domain containing 6B

Gene ID:

51633

## Detection Antibody Information

Catalog Number:  
68450-3-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G Magarose purification

Clone ID:  
3A12A11

Reactivity:  
human

GenBank:  
BC029760

Immunogen Catalog Number:  
Ag33607

Conjugate:

Unconjugated

Full name:

OTU domain containing 6B

Gene ID:

51633

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.098-25 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50962-2 targets OTUD6B in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: OTUD6B Monoclonal antibody, PBS Only (Capture) 68450-1-PBS (1D1G1). 100  $\mu$ g. Concentration 1 mg/mL.

Detection antibody: OTUD6B Monoclonal antibody, PBS Only (Detector) 68450-3-PBS (3A12A11). 100  $\mu$ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50962-2, OTUD6B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68450-1-PBS. Detection antibody: 68450-3-PBS. Standard: Ag33607. Range: 0.098-25 ng/mL