

For Research Use Only

NUS1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50965-3

Capture Antibody Information

Catalog Number:
60672-4-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
2A7B10

Reactivity:
human

GenBank:
BC013026

Immunogen Catalog Number:
Ag8959

Conjugate:
Unconjugated

Full name:
nuclear undecaprenyl pyrophosphate
synthase 1 homolog (*S. cerevisiae*)

Gene ID:
116150

Detection Antibody Information

Catalog Number:
60672-5-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
2F1B2

Reactivity:
human

GenBank:
BC013026

Immunogen Catalog Number:
Ag8959

Conjugate:
Unconjugated

Full name:
nuclear undecaprenyl pyrophosphate
synthase 1 homolog (*S. cerevisiae*)

Gene ID:
116150

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-100 ng/mL (Cytometric Bead
Array)

Recommended Dilutions:
It is recommended that this reagent
should be titrated in each testing
system to obtain optimal results.

Product Information

MP50965-3 targets NUS1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NUS1 Monoclonal antibody, PBS Only (Capture) 60672-4-PBS (2A7B10). 100 μ g. Concentration 1 mg/mL.

Detection antibody: NUS1 Monoclonal antibody, PBS Only (Detector) 60672-5-PBS (2F1B2). 100 μ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50965-3, NUS1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60672-4-PBS. Detection antibody: 60672-5-PBS. Standard:Ag8959. Range: 0.098-100 ng/mL.