

For Research Use Only

# SPAG8 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50977-1

## Capture Antibody Information

<b>Catalog Number:</b> 60679-1-PBS	<b>Clone ID:</b> 1D12E2	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Mouse	<b>Reactivity:</b> human	<b>Full name:</b> sperm associated antigen 8
<b>Isotype:</b> IgG1	<b>GenBank:</b> BC019247	<b>Gene ID:</b> 26206
<b>Purification Method:</b> Protein G Magarose purification	<b>Immunogen Catalog Number:</b> Ag5150	

## Detection Antibody Information

<b>Catalog Number:</b> 60679-2-PBS	<b>Clone ID:</b> 1C6H5	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Mouse	<b>Reactivity:</b> human	<b>Full name:</b> sperm associated antigen 8
<b>Isotype:</b> IgG1	<b>GenBank:</b> BC019247	<b>Gene ID:</b> 26206
<b>Purification Method:</b> Protein G Magarose purification	<b>Immunogen Catalog Number:</b> Ag5150	

## Applications

<b>Tested Applications:</b> Cytometric bead array	<b>Range:</b> 1.563-100 ng/mL (Cytometric Bead Array)	<b>Recommended Dilutions:</b> It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

## Product Information

MP50977-1 targets SPAG8 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SPAG8 Monoclonal antibody, PBS Only (Capture) 60679-1-PBS (1D12E2). 100 µg. Concentration 1 mg/ml.

Detection antibody: SPAG8 Monoclonal antibody, PBS Only (Detector) 60679-2-PBS (1C6H5). 100 µg. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50977-1, SPAG8 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60679-1-PBS. Detection antibody: 60679-2-PBS. Standard: Ag5150. Range: 1.563-100 ng/mL.