

For Research Use Only

FAM38A Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50981-1

Capture Antibody Information

Catalog Number:
60682-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
1A8A5
Reactivity:
human
GenBank:
BC008073
Immunogen Catalog Number:
Ag7791

Conjugate:
Unconjugated
Full name:
family with sequence similarity 38, member A
Gene ID:
9780

Detection Antibody Information

Catalog Number:
60683-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
2E5F1
Reactivity:
human
GenBank:
BC008073
Immunogen Catalog Number:
Ag29858

Conjugate:
Unconjugated
Full name:
family with sequence similarity 38, member A
Gene ID:
9780

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50981-1 targets FAM38A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: FAM38A Monoclonal antibody, PBS Only (Capture) 60682-1-PBS (1A8A5). 100 µg. Concentration 1 mg/ml.
Detection antibody: FAM38A Monoclonal antibody, PBS Only (Detector) 60683-1-PBS (2E5F1). 100 µg. Concentration 1 mg/ml.
Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.
Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50981-1, FAM38A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60682-1-PBS. Detection antibody: 60683-1-PBS. Standard:Ag7887. Range: 0.098-100 ng/mL.