

For Research Use Only

EXT2 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50989-3

Capture Antibody Information

Catalog Number: 60686-4-PBS	Clone ID: 1G6F9	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: exostoses (multiple) 2
Isotype: IgG1	GenBank: BC010058	Gene ID: 2132
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag1748	

Detection Antibody Information

Catalog Number: 60686-2-PBS	Clone ID: 1D1E11	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: exostoses (multiple) 2
Isotype: IgG1	GenBank: BC010058	Gene ID: 2132
Purification Method: Protein G purification	Immunogen Catalog Number: Ag1748	

Applications

Tested Applications: Cytometric bead array	Range: 3.125-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
------------------------------------------------------	----------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------

Product Information

MP50989-3 targets EXT2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: EXT2 Monoclonal antibody, PBS Only (Capture) 60686-4-PBS (1G6F9). 100 µg. Concentration 1 mg/ml.

Detection antibody: EXT2 Monoclonal antibody, PBS Only (Detector) 60686-2-PBS (1D1E11). 100 µg. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50989-3, EXT2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60686-4-PBS. Detection antibody: 60686-2-PBS. Standard: Ag1748. Range: 3.125-100 ng/mL.