

For Research Use Only

AXL Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50990-2

Capture Antibody Information

Catalog Number: 68108-2-PBS	Clone ID: 1A4D11	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: AXL receptor tyrosine kinase
Isotype: IgG1	GenBank: BC032229	Gene ID: 558
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag14401	

Detection Antibody Information

Catalog Number: 68108-4-PBS	Clone ID: 1E10D3	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: AXL receptor tyrosine kinase
Isotype: IgG1	GenBank: BC032229	Gene ID: 558
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag14401	

Applications

Tested Applications: Cytometric bead array	Range: 0.195-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP50990-2 targets AXL in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: AXL Monoclonal antibody, PBS Only (Capture) 68108-2-PBS (1A4D11). 100 µg. Concentration 1 mg/ml.

Detection antibody: AXL Monoclonal antibody, PBS Only (Detector) 68108-4-PBS (1E10D3). 100 µg. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50990-2, AXI Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68108-2-PBS. Detection antibody: 68108-4-PBS. Standard: Ag14401. Range: 0.195-100 ng/mL.